

Perforin Antibody : B-D48 Azide free

ANTI HUMAN PERFORIN AZIDE FREE

CATALOGUE N° : 854.950.000

CLONE : B-D48

ISOTYPE : Mouse IgG1

IMMUNISATION : Natural human perforin

HYBRIDOMA : Myeloma X63/AG.8653 x Balb/c spleen cells

SPECIFICITY : Recognizes the human perforin, a 70 kDa protein

QUANTITY : 200 µg / 200 µl

PURIFICATION : Ion exchange chromatography

FORMAT : Phosphate-buffered saline. Sterile-filtered through 0.22 µm and treated to remove endotoxins. Carrier and preservative free

STORAGE : Stable at +2-8°C for 12 months
For longer storage, freeze aliquots

APPLICATION : Suitable for flow cytometry
Suitable for Western-blot (1 µg/ml)
Suitable for immunohistochemistry on frozen sections (20 µg/ml)

Version 2 – 10.07

For research use only