

## Perforin Antibody : B-D48 FITC conjugated

ANTI HUMAN PERFORIN FITC CONJUGATED

**CATALOGUE N° :** 854.951.010

**CLONE :** B-D48

**ISOTYPE :** Mouse IgG1

**IMMUNISATION :** Natural human perforin

**HYBRIDOMA :** Myeloma X63/AG.8653 x Balb/c spleen cells

**SPECIFICITY :** Recognizes the human perforin, a 70 kDa protein

**QUANTITY :** 100 tests / 1ml

**PURIFICATION :** Ion exchange chromatography

**FORMAT :** Phosphate-buffered saline with 5% BSA and 0.1% sodium azide

**WORKING DILUTION :** Use 10 µl to label 10<sup>6</sup> cells or 100 µl of whole blood

**STORAGE :** Stable at +2-8°C for 12 months.  
For longer storage, freeze aliquots

**APPLICATION :** Suitable for flow cytometry