

## Perforin Antibody : B-D48 Unconjugate

ANTI HUMAN PERFORIN UNCONJUGATE

**CATALOGUE N° :** 854.953.020

**CLONE :** B-D48

**ISOTYPE :** Mouse IgG1

**IMMUNISATION :** Natural human perforin

**HYBRIDOMA :** Myeloma X63/AG.8653 x Balb/c spleen cells

**SPECIFICITY :** Recognizes the human perforin, a 70 kDa protein

**QUANTITY :** 200 tests / 2 ml

**PURIFICATION :** Ion exchange chromatography

**FORMAT :** Phosphate-buffered saline with 1% BSA and 0.1% sodium azide


**WORKING DILUTION :** Use 10 µl to label 10<sup>6</sup> cells or 100 µl of whole blood

**STORAGE :** Stable at +2-8°C for 12 months

**APPLICATION :** Suitable for flow cytometry



*Flow cytometric analysis of perforin in human YT-Indy cell line with the B-D48 monoclonal antibody*

 Sodium azide is a poisonous and hazardous substance, handle and dispose of properly